

GelVue Dual UV Transilluminator-Quick Guide

■ Description

1. GelVue Dual UV Transilluminator is an ultra violet light box with two UV wavelengths, for use in gel imaging.

Model


| Model | Filter Size | Wavelength (nm) |
|-----------|-------------|-----------------|
| GVLM20 | 20x20 | 365 and 302 |
| GVLM30 | 20X30 | 365 and 302 |
| GVLM25/30 | 25x30 | 365 and 254 |
| GVLS20 | 20x20 | 365 and 254 |
| GVLS30 | 20x30 | 365 and 254 |
| GVLS25/30 | 25x30 | 365 and 254 |
| GVSM20 | 20x20 | 254 and 302 |
| GVSM30 | 20x30 | 254 and 302 |
| GVSM25/30 | 25x30 | 254 and 302 |

Key: L- Long (365nm), M- Medium (302nm) and S- short (254nm)

WARNING: UV radiation can be dangerous for unprotected eyes and skin, therefore, appropriate precautions must be taken. UV protective goggles or a face shield must be worn. Please contact our sales office for further details and accessories.

■ Operating instructions

1. Select the wavelength of choice.
2. Make sure that the On/Off switch is in the Off position.



3. Plug the mains input to the unit and supply. DO NOT switch the unit on before checking that safety precautions have been observed.
4. Now switch the unit on and ensure the switch lamp illuminates and the lamps in the unit illuminate evenly. N.B. - only half the lamps should illuminate (every other one).
5. Please switch off the unit when not in use.
6. If unit fails to light up, check that the mains supply is on, that the On/Off switch is in the On position. If no fault can be found please contact our service department for assistance.

CONTACT SYNGENE:
EUROPE:
Tel: +44 (0)1223 727123
Email: support@syngene.com

USA:
Tel: 800 686 4407/301 662 2863
Email: support@syngene.com

www.syngene.com
GVDT0709