

## What is background subtraction?

### Introduction

Baseline correction is used to reduce the effect of background and to emphasise the signal. All automatic background subtraction can be carried out in the integration parameters box, which appears when you click on the icon found on the GeneTools header bar:



**Integration parameters icon**

The automatic background options include:

- None
- Track Borders
- Lowest Slope
- Track Borders & Slope
- Rolling Disk

### None

No background corrections are made.

### Track Borders

The signal is calculated relative to the points at the track border. This means that the first pixel outside the track is used for background correction of the other pixels inside. The outside value is probably an average of the first pixels either side of the track.

### Lowest Slope


To understand how lowest slope works consider the following example, which goes through the process as if you were doing it by hand. GeneTools automatically calculates this for you and corrects the background (**Figure 1e**).

To find the first segment of the baseline, the program starts at one end of the profile and draws imaginary lines from the start point to all other points on the profile (**Figure 1b**).

The program then chooses the line of lowest slope. This is the red line (**Figure 1c**) as the green and blue lines have larger (less negative) slopes (lines to all other points on the profile have even larger slopes). To find the second segment of the baseline, the program moves to the end of the first segment and again draws imaginary lines from here to the remaining points on the profile.

The software again chooses the line of lowest slope (the red line) and repeats this process until it gets to the other end of the profile (**Figure 1d**).

The corrected signal is then measured as the height above this baseline (**Figure 1e**)



**Figure 1- Example of lowest slope background subtraction**  
In this example only three lines are shown for simplicity.

### Track Borders & Slope

A combination of Track borders and lowest slope whereby the signal is corrected for the track borders, and then the lowest slope correction is made.

### Rolling Disk

For this method the program first calculates the position of the line formed by the centre of a disk with the set radius, rolled along below the profile.

The baseline is then one radius length above this line, and the corrected signal is measured as the height above this baseline. The larger the radius of the disc, the less the background rises with the profile.




Figure 2- Background subtraction using the rolling disk method

### Manual Baseline Corrections

You can use the Integration parameters dialog box to set the type of automatic baseline correction to be used for a sample, but you can also set a manual baseline.

To set a manual baseline for a track:

1. Click on the track in the Image pane to select it.
2. Choose Edit manual baseline from the Track menu or the context menu displayed when you right-click in the Image pane or Profile pane. The baseline will be displayed with circular drag handles at its ends and in the colour selected for the manual baseline (Figure 3a)
3. If required, double-click on the baseline to add additional drag handles so that you can bend the line. **N.B.** Any peak markers and boundaries will be removed from the profile (Figure 3b).
4. To re-shape the baseline, position the pointer over the drag handle you want to move. The pointer will change to a four-way arrow (Figure 3c). Drag the handle to reshape the baseline (Figure 3d) then drop the handle in its new position (Figure 3e).

**N.B.** While you are dragging a handle, the new baseline may be drawn above the profile, but when you release the mouse button it will be reshaped so that it always lies below it.




Figure 3- Setting a manual baseline for a track

**Syngene reserves the right to amend or change specifications without prior notice. This Application Note supersedes all earlier versions.**

All trademarks acknowledged.

May 2010

UK tel: +44 (0)1223 727123  
Email: [sales@syngene.com](mailto:sales@syngene.com)

USA tel: 800 686 4407/301 662 2863  
Email: [ussales@syngene.com](mailto:ussales@syngene.com)