

U:Genius³ EZ Installation - Quick Guide

■ System components

1. Darkroom (Epi White lighting factory fitted)
2. 3MP camera with zoom lens pre-fitted to camera
3. Stylus which can be attached to the darkroom in a suitable location
4. Mains power cable and 1 metre USB cable

CAUTION: Do not connect the mains power cable to the UG3 until you are satisfied that everything is connected correctly. For assistance please contact your supplier or Syngene directly


■ Installation

Transilluminator set-up

1. Place the UV transilluminator inside and connect it to the free flying mains lead inside the U:Genius³ EZ system

Camera set-up

1. Connect the USB cable from the camera and the serial cable from the zoom-lens to the back of the U:Genius³ EZ unit



Installing the filter

1. Place the UV filter with the thread pointing upwards in the filter drawer (found below the camera and on the top of the U:Genius³ EZ system)
2. This filter can be interchanged with other filters for different dyes



Power on

1. Connect the mains lead from the U:Genius³ EZ system to the mains power supply and switch on
2. Power on by **momentarily** pressing the round black button on the rear of the unit. There will be a 1-2 second delay before a picture appears on the display



■ Select Regional Power Frequency

1. From the first screen select 'Preferences' icon
2. Select either 50 or 60Hz



or



■ Installing GeneTools on a separate PC

1. Plug the Syngene branded flash drive into the PC you wish to install GeneTools on. Navigate into the "GeneTools" folder on the flash drive and run the "InstallGeneTools" program
2. Follow the on-screen instructions to install and license the GeneTools software. You will probably have been provided with media keys to license the software. These media keys can be found in a file named "Media Keys.txt" on the Syngene branded flash drive

CONTACT SYNGENE:
EUROPE:
Tel: +44 1223 727123
Email: support@syngene.com

USA:
Tel: 800 686 4407 or +1 301 662 2863
Email: support@syngene.com

www.syngene.com
UG3EZINST0416