

Recommended environment for Syngene image capture systems

The recommended operating specifications are as follows:

1. **Operating temperature:** 5 to 35 degrees centigrade
2. **Operating humidity:** 15 to 80% at 26 degrees centigrade (non-condensing)
3. **Storage temperature:** -10 to 65 degrees centigrade
4. **Storage humidity:** 90% at 65 degrees centigrade (non-condensing)

Suggested ambient conditions are 21 degrees centigrade with a relative humidity of 45% as this is an easy level to maintain.

This level provides the best protection from the static discharge common at lower humidity levels and the corrosive effects of higher humidity levels. These conditions also provide good operator comfort levels.

Syngene reserves the right to amend or change specifications without prior notice. This Application Note supersedes all earlier versions.

All trademarks acknowledged.

31.08.09

UK tel: +44 (0)1223 727123
Email: sales@syngene.com

USA tel: 800 686 4407/301 662 2863
Email: ussales@syngene.com