



## DigiGenius

| Part number  | Description   | Notes   |
|--|---|---|
| <b>DigiGenius System</b>   | With hood, UV filter & 1 copy of GeneTools. Transilluminator & converter not supplied please select required source from list below |   |
| DG1  | DigiGenius system, with 10M Pixel camera  | Not voltage dependant   |
| <b>Stand-alone Single wavelength slide-out transilluminators</b> |   | For use with DG or on Benchtop as stand-alone. Includes mains cable |
| BT-2020L-E   | Transilluminator (20 x 20cm; 365nm) 230V, 50Hz  | Longwave (no cover) 6 tubes   |
| BT-2020L-E60   | Transilluminator (20 x 20cm; 365nm) 230V, 60Hz  | Longwave (no cover) 6 tubes   |
| BT-2020L-A   | Transilluminator (20 x 20cm; 365nm) 115V, 60Hz  | Longwave (no cover) 6 tubes   |
| BT-2020L-J   | Transilluminator (20 x 20cm; 365nm) 100V, 50/60Hz   | Longwave (no cover) 6 tubes   |
| BT-2020M-E   | Transilluminator (20 x 20cm; 302nm) 230V, 50Hz  | Mediumwave (no cover) 6 tubes                                       |
| BT-2020M-E60   | Transilluminator (20 x 20cm; 302nm) 230V, 60Hz  | Mediumwave (no cover) 6 tubes                                       |
| BT-2020M-A   | Transilluminator (20 x 20cm; 302nm) 115V, 60Hz  | Mediumwave (no cover) 6 tubes                                       |
| BT-2020M-J   | Transilluminator (20 x 20cm; 302nm) 100V, 50/60Hz   | Mediumwave (no cover) 6 tubes                                       |
| BT-2020S-E   | Transilluminator (20 x 20cm; 254nm) 230V, 50Hz  | Shortwave (no cover) 6 tubes  |
| BT-2020S-E60   | Transilluminator (20 x 20cm; 254nm) 230V, 60Hz  | Shortwave (no cover) 6 tubes  |
| BT-2020S-A   | Transilluminator (20 x 20cm; 254nm) 115V, 60Hz  | Shortwave (no cover) 6 tubes  |
| BT-2020S-J   | Transilluminator (20 x 20cm; 254nm) 100V, 50/60Hz   | Shortwave (no cover) 6 tubes  |

|                       |  |   |
|-----------------------|--|---|
|                       |  |   |
| <b>Filters</b>        |  |   |
| FILTUV-58             | UV filter (UV032)  | For IG/UG & DG units Included                         |
| FILTSP-58             | Short pass filter (SW032)                                  | For IG/UG & DG units                                  |
| FILTLP-58             | Long pass filter (LW032)                                   | For IG/UG & DG units                                  |
|                       |  |   |
| <b>Options</b>        |  |   |
| GE-CONVERT5           | Visible light converter                                    | For IG/UG & DG systems                                |
|                       |  |   |
| ULTRA-SLIM-LED        | Blue slimline LED transilluminator, 10x12cm + sliding tray | For IG/UG & DG systems. Not voltage dependant         |
| GENEDIRECTORY-DB      | GeneDirectory Database software                            |   |
| P93DW                 | Mitsubishi Digital Thermal Printer                         |   |
| K65HM                 | Thermal paper, Matt, 1 roll                                | For thermal printers                                  |
| K91HG                 | Thermal paper, Glossy, 1 roll                              | For thermal printers                                  |
| UV-SCREEN             | Safety screen  | For standing in front of transilluminator or darkroom |
|                       |  |   |
| <b>Spare UV Tubes</b> |  |   |
| SWT8                  | Short wave tube, 8W, 254nm                                 |   |
| MWT8                  | Medium wave, 8W, 302nm                                     |   |
| LWT8                  | Longwave 8W, 365nm   |   |