



UGenius 3

Part number	Description	Notes
U:GENIUS 3 System	With epi LED white lights, filter drawer & UV filter. Integral processor so no PC required. Large 8.7 inch Touch screen monitor for system control. With image capture software & 1 copy of GeneTools. Max image size 20x20cm. Transilluminator & converter not supplied please select required source from list below	
UGEN3	UGenius 3MP 12/16bit system	Not voltage dependant
Slimline Single wavelength slide-out transilluminators		For use with IG and UG systems
GE-2020L-E	Transilluminator (20 x 20cm; 365nm) 230V, 50Hz with runners	Longwave - with runners (no cover) 6 tubes
GE-2020L-E60	Transilluminator (20 x 20cm; 365nm) 230V, 60Hz with runners	Longwave - with runners (no cover) 6 tubes
GE-2020L-A	Transilluminator (20 x 20cm; 365nm) 115V, 60Hz with runners	Longwave - with runners (no cover) 6 tubes
GE-2020L-J	Transilluminator (20 x 20cm; 365nm) 100V, 50/60Hz with runners	Longwave - with runners (no cover) 6 tubes
GE-2020M-E	Transilluminator (20 x 20cm; 302nm) 230V, 50Hz with runners	Mediumwave - with runners (no cover) 6 tubes
GE-2020M-E60	Transilluminator (20 x 20cm; 302nm) 230V, 60Hz with runners	Mediumwave - with runners (no cover) 6 tubes
GE-2020M-A	Transilluminator (20 x 20cm; 302nm) 115V, 60Hz with runners	Mediumwave - with runners (no cover) 6 tubes
GE-2020M-J	Transilluminator (20 x 20cm; 302nm) 100V, 50/60Hz with runners	Mediumwave - with runners (no cover) 6 tubes
GE-2020S-E	Transilluminator (20 x 20cm; 254nm) 230V, 50Hz with runners	Shortwave - with runners (no cover) 6 tubes
GE-2020S-E60	Transilluminator (20 x 20cm; 254nm) 230V, 60Hz with runners	Shortwave - with runners (no cover) 6 tubes
GE-2020S-A	Transilluminator (20 x 20cm; 254nm) 115V, 60Hz with runners	Shortwave - with runners (no cover) 6 tubes
GE-2020S-J	Transilluminator (20 x 20cm; 254nm) 100V, 50/60Hz with runners	Shortwave - with runners (no cover) 6 tubes

Filters 58mm diameter		
FILTUV-58	UV filter (UV032)	For IG/UG & DG units Included
FILTSP-58	Short pass filter (SW032)	For IG/UG & DG units
FILTLP-58	Long pass filter (LW032)	For IG/UG & DG units
Options		
GE-CONVERT5	Visible light converter	For IG/UG & DG systems
ULTRA-SLIM-LED	Blue slimline LED transilluminator, 10x12cm + sliding tray	For IG/UG & DG systems. Not voltage dependant
GE-NFS	Neutral field screens (1 blue, 1 frosted).	Corrects for any uneven Illumination by acquiring a "neutral-field" reference image with the same intensity illumination as the experiment. For use with Epi-UV lighting and Transilluminators in IG & UG systems
MPVIEWER	Microplate Viewer Anti parallax tray	For transillumination of microtitre plates. Not for use with epi lighting
GENEDIRECTORY-DB	GeneDirectory Database software	
P93DW	Mitsubishi Digital Thermal Printer	
K65HM	Thermal paper, Matt, 1 roll	For thermal printers
K91HG	Thermal paper, Glossy, 1 roll	For thermal printers
UV-SCREEN	Safety screen	For standing in front of transilluminator or darkroom
Spare UV Tubes		
SWT8	Short wave tube, 8W, 254nm	
MWT8	Medium wave, 8W, 302nm	
LWT8	Longwave 8W, 365nm	